DDO2 Portable Oxygen Analyzer

Instruction Manual

Bedienungsanleitung / Manuel d'instructions /
Manual de instrucciones / Manuale di istruzioni

English	1
German	15
French	29
Spanish	43
Italian	57

Contents

Contents	1
Introduction	2
About this Product	2
About this Manual	2
Precautions for Use	3
Cautions	3
Description of parts and accessories.....	4
Appearance of the main unit	4
Display screen	5
Accessories.....	5
Instructions for use	6
Powering on	6
Low battery	6
Measuring oxygen concentration and flow rate.....	6
Pressure Measurement.....	7
Powering off.....	7
Replacing the batteries	8
Verification.....	9
Maintenance and Storage	10
Maintenance	10
Storage.....	10
Calibration.....	10
Troubleshooting	11
Safety and Regulatory Information	12
Specifications	13
Warranty and Service	14
Product Warranty.....	14
After-sales service	14

Introduction

About this Product

The DDO2 is a dedicated measuring instrument for oxygen concentrators and is an essential piece of equipment for home oxygen therapy services. It allows you to easily measure oxygen concentration, flow rate and outlet pressure all in one unit. Thanks to the use of an ultrasonic sensor, there is no need for on-site atmospheric adjustment or special training; operation is simple and the device is ready for immediate use. The DDO2 has the following features.

Intuitive Operation

The single-button design ensures intuitive operation. Anyone can operate it easily, so even first-time users can take measurements without hesitation. This saves time and allows you to obtain measurement results quickly.

High-Visibility Display

The OLED display enhances visibility and clearly shows measurement results. Even in dimly lit environments during oxygen concentrator inspections, the necessary information can be confirmed at a glance, so no time is wasted reading the results.

No Atmospheric Conditioning Required

No tedious atmospheric conditioning is required. Thanks to the ultrasonic oxygen sensor, measurements can begin immediately after switching on the power. Accurate measurements of oxygen concentration, flow rate and pressure are possible without any waiting time. This contributes to more efficient inspection work.

Cost Reduction

Three functions—oxygen concentration, flow rate and outlet pressure—are consolidated into a single unit. There is no need to manage multiple measuring instruments. The ultrasonic sensor requires no replacement, ensuring low running costs. We handle calibration in-house, and calibration certificates are available as an optional service.

About this Manual

- Please use this device correctly, observing the precautions outlined in this manual.
- Please keep this instruction manual in a safe place so that you can refer to it at any time.

Precautions for Use


Cautions

- This product is intended for use only by medical professionals or those with specialist knowledge, such as home healthcare service providers. It is not intended for use by patients or the general public.
- Do not use this product for any purpose other than with an oxygen concentrator. This product is a measuring instrument designed for use with oxygen concentrators.
- Do not apply pressure exceeding 250 kPa to this product. This may cause the sensor to malfunction, resulting in inaccurate measurements.
- Do not subject the product to strong impacts or vibrations.
- This product will not display correct readings when measuring pure oxygen.
- Do not breathe into the measurement port. This will cause high humidity, preventing the display of correct readings.
- Allow the meter to acclimatise to room temperature before use.
- Insert the battery in the specified direction.
- Remove the batteries if the device is not to be used for a long period of time.
- Please use within the specified operating temperature and humidity range.
- Use near strong magnetic fields or high-frequency equipment may cause measurement errors.
- When cleaning or disinfecting, take care to ensure that no liquid enters the measurement port or electronic components.
- Ensure that water, other liquids, dust or fine particles do not come into contact with or enter the instrument.
- Take care to avoid dropping the device or subjecting it to strong impacts.
- Do not store in locations subject to high temperatures, high humidity or direct sunlight.

Description of parts and accessories

Appearance of the main unit

Display screen

Oxygen concentration
Displayed as in %.

Flow Rate
Displayed in LPM.

Oxygen concentration and flow rate
measurement screen

Low battery icon


When the battery level
is low.

Pressure
Displayed in kPa.

Pressure measurement screen

Accessories

Case cover
A silicone cover Silicone
cover.

Connection Tube
A silicone tube used to
connect this product to the
oxygen concentrator.

Instruction manual
This manual.

Instructions for use

Powering on

- 1) When you press the POWER/MODE button, the Drive Medical logo screen will be displayed for 3 seconds.
- 2) Once the screen displays the oxygen concentration and flow rate, the device is ready for use.

Low battery


If the battery icon flashes as shown in the figure on the right when the power is switched on, this indicates that the battery is low. Please replace the battery with a new one.

Measuring oxygen concentration and flow rate

- 1) Insert the connection tube into the measurement port and connect the other end to the oxygen concentrator to be measured.
- 2) The measured oxygen concentration and flow rate will be displayed simultaneously on the screen.

Pressure Measurement

- 1) Press the POWER/MODE button on the oxygen concentration and flow rate measurement screen to switch to the pressure measurement screen.
- 2) Pressing the outlet with your finger will cause the pressure reading to rise; please do not remove your finger until the reading stabilises.
- 3) Once the pressure reading has stabilised, read the measurement. The reading will return to 0 kPa when you remove your finger from the outlet.

Powering off

To switch off the power, press and hold the POWER/MODE button for 3 seconds. The unit will also switch off automatically if the flow rate display remains at 0.0 (LPM) for 2 minutes.

Replacing the batteries

- 1) Remove the main unit from the case cover.

It can be easily removed by pulling it out from the bottom of the unit.

- 2) Press the ▽ mark on the battery cover on the back of the unit with your finger and slide it in the direction of the arrow to remove it.

- 3) Remove the old batteries and replace them with new ones (two AA batteries), ensuring the polarity matches the markings on the bottom of the battery compartment.

- 4) Slide the battery cover back into place to close it.





- 5) Place the main unit back into the case cover.

Verification

To verify that the device remains calibrated, please carry out a performance check by following the steps below.

- 1) Connect a source of pure oxygen (>99.9% O₂) to the measurement port and set the flow rate to 3 – 5 L/min and the pressure to 200 kPa or less.
- 2) With the power off, press and hold the POWER/MODE button for 5 seconds.
- 3) The verification mode will start and measurement will begin, so start the flow of oxygen.
- 4) If the measured value is within the range of 98.0 % to 102.0 % (100 %O₂ ±2.0 %), the unit is operating normally.

If the reading falls outside this range, "Err" will flash to indicate a fault.

If "Err" continues to flash even after repeating steps 2 to 4, the unit may require repair. Please contact your dealer or distributor.

- Please use oxygen gas with a purity of 99.9% or higher as the source of pure oxygen.
- Performance cannot be verified correctly using gas from an oxygen concentrator.

Maintenance and Storage

Maintenance

- Wipe the surface of the unit regularly with a dry, clean cloth.
- This device is not waterproof. Please take care to ensure it does not come into contact with water.

Storage

- After use, keep the device clean and protect it from moisture and dust. Also, do not leave it in extreme temperatures.
- If you do not intend to use this device for an extended period, please remove the batteries.

Calibration

This product is designed to be calibration-free, and regular calibration is not required under normal conditions of use.

However, performing regular calibration allows you to verify measurement accuracy and assess the condition of the device, ensuring reliable long-term use.

Calibration is carried out by the manufacturer or an organisation certified by the manufacturer; if you wish to have the device calibrated, please contact the retailer from whom you purchased it.

Troubleshooting

If you suspect a fault, please check the following before requesting an inspection or repair. If the problem persists, please stop using the product and contact your supplier to arrange for an inspection or repair.

Symptoms	Cause	Solution
The power does not turn on	Low battery	Please replace the batteries.
	The battery is inserted the wrong way round	Please insert the batteries with the correct polarity.
Oxygen concentration is not displayed correctly	Not correctly connected to the oxygen concentrator	Please check that it is securely connected to the oxygen concentrator.
	The connecting tube is cracked or damaged	Check the tubing for cracks or damage and replace it if any abnormalities are found.
Flow rate cannot be measured correctly	Not correctly connected to the oxygen concentrator	Please ensure that it is securely connected to the oxygen concentrator.
	The connecting tube is cracked or damaged	Check the tube for cracks or damage, and replace it if any abnormalities are found.
Pressure cannot be measured correctly	The measurement method is incorrect	When taking a measurement, ensure that the outlet is securely covered with your finger.
	Not correctly connected to the oxygen concentrator	Please check that it is securely connected to the oxygen concentrator.
	The connecting tube is cracked or damaged	Check the tubing for cracks or damage, and replace it if any abnormalities are found.
Oxygen concentration and flow rate readings are displayed as "--"	The flow rate exceeds the measurement range	Please adjust the oxygen flow rate to within the measurement range. The flow rate measurement range for this product is a maximum of 10 LPM.

Safety and Regulatory Information

Manufacturer Information

Nihon Rufuto Co., Ltd.

7-2-1 Uchiya, Minami-ku, Saitama-shi, Saitama 336-0034, Japan

www.nihon-rufuto.com

Authorised Representative / Importer within the EU

Drive Medical GmbH & Co. KG

Leutkircher Str. 44 | 88316 Isny | Germany

www.drivedevilbiss.com

Compliance with the EMC Directive and RoHS Directive

This product complies with the European Union (EU) Electromagnetic Compatibility Directive (EMC Directive 2014/30/EU) and the Directive on the Restriction of the Use of Certain Hazardous Substances (RoHS Directive 2011/65/EU).

Explanation of symbols

	Caution		Minor caution, information
	Product model number		Product serial number
	Manufacturer of this product		Authorised Representative within the EU
	This indicates that this product complies with the applicable EU directives. This indicates that this product complies with the applicable EU directives.		
	Do not dispose of this product as general waste. Dispose of it in accordance with the WEEE Directive.		
	TRIMAN MARK		

Specifications

Product Name	DDO2 Portable Oxygen Analyzer	
Model	HC2-A-101	
Dimensions	76 x 35 x 135 mm (case cover, excluding protrusions)	
Weight	236 g (excluding case cover)	
Detection method	Ultrasonic	
Oxygen concentration	Measurement range	21.0 – 95.6 %
	Accuracy	±2.0 %FS
	Resolution	0.1 %
Flow rate	Measuring range	0 – 10 L/min
	Accuracy	±0.2 L/min
	Resolution	0.1 L/min
Pressure	Measuring range	0 – 200 kPa
	Accuracy	±2 kPa
	Resolution	1 kPa
Response time	1 second	
Start-up time	3 seconds	
Operating environment	Temperature	5 – 35 °C
	Humidity	30 – 85 %RH (no condensation)
Storage environment	Temperature	-20 – 50 °C
	Humidity	0 – 95 %RH (no condensation)
Auto-off function	2 minutes	
Power supply	2 x AA batteries	
Alarm Indication	Low battery	
Accessories	Case cover, Connecting tube, Instruction manual	

Warranty and Service

Product Warranty

The warranty period is one year from the date of purchase. Should a fault occur with the device when used in accordance with the instructions in the user manual, we will repair it free of charge (including parts replacement and labour costs). However, repairs will be charged in the following cases:

- If a fault occurs due to use not in accordance with the manual
- Damage caused by wilful misconduct or negligence

Adjustments and repairs will be carried out by our Technical Department. Please note that for adjustments and repairs after the warranty period, should the cause be confirmed, a chargeable repair will be carried out following the submission of an estimated quotation. Although the equipment is calibrated for performance at the time of dispatch, please contact the enquiry details below should you experience any abnormalities, faults or damage.

After-sales service

Please contact the retailer from whom you purchased the product regarding calibration and repairs.

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