




## Portable Oxygen Concentrator

Model: OX-2A

User Manual



# TABLE OF CONTENTS

<b>1. Foreword</b> .....	<b>4</b>
<b>1.1 General Information</b> .....	<b>4</b>
<b>1.2 Classification</b> .....	<b>4</b>
<b>1.3 Typographical Conventions</b> .....	<b>4</b>
<b>2. Principle of operations and Intended Use</b> .....	<b>5</b>
<b>2.1 Principle of Operations</b> .....	<b>5</b>
<b>2.2 Intended Use</b> .....	<b>5</b>
<b>3. Safety Instruction</b> .....	<b>7</b>
<b>3.1 Intended user profile</b> .....	<b>7</b>
<b>3.2 Warnings Overview</b> .....	<b>7</b>
<b>3.3 Cautions Overview</b> .....	<b>10</b>
<b>3.4 Overview of Important Information</b> .....	<b>11</b>
<b>4. Instructions and Training</b> .....	<b>11</b>
<b>5. Product Description</b> .....	<b>13</b>
<b>5.1 Device and Accessories Description</b> .....	<b>13</b>
<b>6. General Instructions Before Use</b> .....	<b>14</b>
<b>6.1 Accessories List</b> .....	<b>14</b>
<b>6.2 Battery</b> .....	<b>14</b>
<b>6.3 Nasal Cannula</b> .....	<b>15</b>
<b>6.4 Humidifier</b> .....	<b>15</b>
<b>6.5 Pull Cart</b> .....	<b>16</b>
<b>7. Operating OX-2A</b> .....	<b>17</b>
<b>7.1 Before Operating Your Oxygen Concentrator</b> .....	<b>17</b>
<b>7.2 Connecting Nasal Cannula</b> .....	<b>18</b>
<b>7.3 Device Overview</b> .....	<b>19</b>
<b>7.4 Control Panel</b> .....	<b>20</b>
<b>7.5 Run Time Screen Description</b> .....	<b>20</b>
<b>7.6 Turning On</b> .....	<b>21</b>
<b>7.7 Adjusting Setting</b> .....	<b>21</b>
<b>7.8 Responding to Alarms</b> .....	<b>22</b>
<b>7.9 Turning Off</b> .....	<b>23</b>
<b>8. Alarm</b> .....	<b>24</b>
<b>9. Troubleshooting</b> .....	<b>27</b>
<b>10. Maintenance and Cleaning</b> .....	<b>28</b>
<b>10.1 Routine Maintenance</b> .....	<b>28</b>
<b>10.2 Cleaning</b> .....	<b>28</b>
<b>10.3 Service Life</b> .....	<b>29</b>
<b>11. Device Repair and Disposal</b> .....	<b>29</b>
<b>11.1 Repair</b> .....	<b>29</b>
<b>11.2 Disposal</b> .....	<b>29</b>
<b>12. Warranty</b> .....	<b>29</b>
<b>13. Disclaimer</b> .....	<b>29</b>
<b>13.1 Disclaimer</b> .....	<b>29</b>
<b>13.2 This Document</b> .....	<b>29</b>
<b>14. Specifications</b> .....	<b>31</b>
<b>15. Traveling With Your OX-2A</b> .....	<b>33</b>
<b>15.1 By Motor Vehicle</b> .....	<b>33</b>
<b>15.2 By Bus or Train</b> .....	<b>33</b>
<b>16. Symbols</b> .....	<b>34</b>
<b>17. EMC Information</b> .....	<b>36</b>

**18. Declaration of Conformity.....38**  
**19. Reporting System of Serious Events..... 39**

## 1. Foreword

Please refer to this manual for detailed instructions on warnings, cautions, specifications, and additional information.

**⚠ IMPORTANT:** Users should read this entire manual before operating the OTM Portable Oxygen Concentrator. Failure to do so could result in personal injury and/or death. If you have questions about the information in this user manual or about the safe operation of this system, contact your distributor.

### 1.1 General Information

This user manual provides information for users of the OXTM Portable Oxygen Concentrator. For the sake of brevity, the terms “concentrator,” “POC,” “unit,” or “device” are sometimes used in this document to refer to the OXTM Portable Oxygen Concentrator. “Patient” and “User” are used interchangeably.

### 1.2 Classification

This device is listed with an internationally recognized testing laboratory and classified with respect to electric shock, fire, and mechanical hazards in accordance with the following standards:

- EN 60601- 1:2006+A2:2021, Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance.
- EN 60601- 1- 2:2015/A1:2021, Part 1- 2: General Requirements for Basic Safety and Essential Performance – Collateral Standard: Electromagnetic Compatibility – Requirements and Tests.
- EN 60601- 1-6:2010+A1:2015 Medical Electrical Equipment – Part 1-6: General Requirements for Basic Safety and Essential Performance – Collateral Standard: Usability.
- EN 60601- 1-8:2007+A1:2013 Medical Electrical Equipment – Part 1-8: General Requirements for Basic Safety and Essential performance – Collateral Standard: General Requirements, Tests and Guidance for Alarm Systems in Medical Electrical Equipment and Medical Electrical Systems.
- EN 60601- 1- 11:2015 General Requirements for Basic Safety and Essential Performance - Collateral Standard: Requirements for Medical Electrical Equipment and Medical Electrical Systems Used in the Home Healthcare Environment.
- Medical Device Regulation (EU) 2017/745.

#### **This equipment is classified as:**

- Class II
- Class IIb according to the REGULATION (EU) 2017/745
- Type BF
- IP22 with the carry bag

### 1.3 Typographical Conventions

This user manual contains warnings, cautions, and notes to help call attention to the most important safety and operational aspects of the device. To help identify these items when they occur in the text, they are shown using the following typographical conventions:

**⚠ WARNING:** Statements that describe serious adverse reactions and potential safety hazards.

**⚠ CAUTION:** Statements that call attention to information regarding any special care to be exercised by the practitioner and/ or patient for the safe and effective use of the device.

**⚠ IMPORTANT:** Statements calling attention to additional significant information about the device or a procedure.

## 2. Principle of Operations and Intended Use

### 2.1 Principle of Operations

The Oxygen Concentrator is a psa oxygen generator. The psa (pressure swing adsorption) oxygen generator is an air separation unit which is based on molecular sieves selective adsorption of nitrogen from air. It can continuously generate oxygen with the purity 90%±3% at ambient temperature. The adsorbed nitrogen can be desorbed by decreasing the adsorbed bed pressure resulting cyclical adsorption-desorption operation.

Pressure Swing Adsorption (PSA) processes rely on the fact that under pressure gases tend to be attracted to solid surfaces, or adsorbed. The higher the pressure, the more gas is adsorbed; when the pressure is reduced, the gas is released, or desorbed. PSA processes can be used to separate gases in a mixture because different gases tend to be attracted to different solid surfaces more or less strongly. If a gas mixture such as air, for example, is passed under pressure through a vessel containing an adsorbent bed that attracts nitrogen more strongly than it does oxygen, part or all of the nitrogen will stay in the bed, and the gas coming out of the vessel will be enriched in oxygen. When the bed reaches the end of its capacity to adsorb nitrogen, it can be regenerated by reducing the pressure, thereby releasing the adsorbed nitrogen. It is then ready for another cycle of producing oxygen enriched air.

### 2.2 Intended Purpose

OXTM portable oxygen concentrator is intended to provide supplemental low flow oxygen for patients suffering from COPD.

The device is portable, enabling patients who need an oxygen device to be treated at home according to a clinician's prescription or direction.

OXTM is not intended for use in life supporting or life sustaining situations, and is provided non-sterile. It is a prescription only device, and it is designed for indoor use, as well as in motor vehicles, buses or train.

### 2.3 Intended Use

OXTM portable oxygen concentrator is intended to provide supplemental low flow oxygen for patients suffering from COPD.

The device is portable, enabling patients who need an oxygen device to be treated at home according to a clinician's prescription or direction.

OXTM is not intended for use in life supporting or life sustaining situations, and is provided non-sterile. It is a prescription only device, and it is designed for indoor use, as well as in motor vehicles, buses or train.

#### **Contraindications:**

OXTM Portable Oxygen Concentrator is not intended to be used:

- in life-supporting or life-sustaining situations
- in an operating or surgical environment
- with a non-adult population
- in conjunction with flammable anaesthetic or flammable materials

#### **Side-effect:**

- Hypercapnia
- Pulmonary toxicity

#### **Intended Patient's specification:**

- a) Population: adult
- b) Age: above 18 years old
- c) Health status: suffering from stage III (severe) or IV (very severe) of COPD with hypoxemia at rest or without hypoxemia.

**Intended target groups:**

Patients suffer from stage III (severe) or IV (very severe) of COPD with hypoxemia at rest or without hypoxemia.

**Medical condition:**

Medical fields concerned is providing supplemental low flow oxygen for patients suffering from at stage III (severe) or IV (very severe) of COPD with hypoxemia at rest or without hypoxemia.

**Clinical benefit:**

The portable oxygen concentrator is intended to provide supplemental low flow oxygen therapy for patients suffering from COPD and correct use of Portable oxygen concentrator can benefit COPD patients by improving their exercise ability to walk more than 70m and by mitigating dyspnoea.

### 3. Safety Instructions

#### 3.1 Intended user profile

**Target user group:** patients and care givers who are trained by an experienced person who has been authorized by the manufacturer and has appropriate training, knowledge, and experience.

Operator	
Age	-Adult (above 18 years old)
Knowledge	minimum: -Read and understand text and Arabic numerals; -Read this manual.
Linguistic	-English or local official language
Education	-At least 18 years old and 12 years intensive reading experience (school). -No maximum.
Experience	- be trained by an experienced person who has been authorized by the manufacturer - And has appropriate training, knowledge and experience.
Permissible impairments	-Mild reading vision impairment or vision corrected to log MAR 0,2(6/10 or 20/32). -Impaired by 40% resulting in 60% of normal hearing at 50 Hz to 2 kHz.

#### 3.2 Warnings Overview

- 1 The device must be used in the carry bag to provide protection from liquid intrusion from rain and/or spills.
- 2 OX-2A for single patient use only.
- 3 There is a risk of fire associated with oxygen equipment and therapy. Do not use near sparks or open flames.
- 4 The settings of OXTM Portable Oxygen Concentrator OX- 2A might not correspond with continuous flow oxygen.
- 5 The settings of other models or brands of portable oxygen concentrators do not correspond with the settings of OXTM portable Oxygen Concentrator OX- 2A.
- 6 Wind or strong drafts can adversely affect accurate delivery of oxygen therapy.
- 7 Geriatrics or any other patient unable to communicate discomfort can require additional monitoring to avoid harm.
- 8 Smoking (including e-cigarettes) during oxygen therapy is dangerous and is likely to result in facial burns, serious injury or death of the patient and others from fire. Do not allow smoking or open flames within the same room as the portable oxygen concentrator or any oxygen carrying accessories. If you smoke, you must always turn the oxygen concentrator off, remove the cannula and leave the room where either the cannula or the concentrator is located. If unable to leave the room, you must wait 10 minutes after the flow of oxygen has been stopped.
- 9 Use only water based lotions that are oxygen compatible, before and during oxygen therapy. Never use petroleum or oil based lotions or salves when operating the device to avoid the risk of fire and burns.
- 10 Open flames during oxygen therapy are dangerous and are likely to result in fire or death. Do not allow open flames within 2 meters of the oxygen concentrator or any oxygen carrying accessory.

- 11 Oxygen makes it easier for a fire to start and spread. Do not leave the nasal cannula on bed coverings or chair cushions with the concentrator on, but not in use; the oxygen will make the materials flammable. Turn the concentrator off when not in use to prevent oxygen enrichment.
- 12 Critical! Explosion hazard. Do not use in the presence of flammable anaesthetics!
- 13 Do not use this device in the presence of pollutants or fumes.
- 14 Do not submerge this device in liquid. Do not expose to water or precipitation. Do not expose to dusty conditions.
- 15 Do not use a device or any accessory that shows any sign of damage.
- 16 Do not use lubricants on this device or any of its accessories.
- 17 Use of this device at an altitude above 2,700 m (9,000 feet), or outside the temperature range of 10 °C (50°F) to 40°C (104°F), or outside the humidity range of 5% to 93% may adversely affect the flowrate and concentrator of oxygen and consequently the quality of therapy. When not in use, the device should be stored in a clean, dry environment between - 20°C and 50°C (-4°F and 120°F). Use and/or storage outside of the valid conditions may damage the product.
- 18 If feeling ill or experiencing discomfort while using this device, contact your clinician or seek medical assistance immediately to avoid harm.
- 19 Your home oxygen provider must verify the compatibility of the device and all accessories used prior to use. To ensure you are receiving the therapeutic amount of oxygen for your medical condition, the device and accessories must only be used after one of more settings have been determined or prescribed for you at your specific activity levels by a healthcare professional.
- 20 The electrical cord and tubing could present a tripping or strangulation hazard. Keep away from children and pets.
- 21 Do not disassemble or modify this device or any of its accessories. Disassembly can create an electric shock hazard and will void the warranty. Contact your distributor for servicing by authorized personnel.
- 22 Use only spare parts recommended by the manufacturer to ensure proper function and to avoid the risk of fire and burns.
- 23 The operator should read and understand this entire manual before using the device.
- 24 The device is not intended for life support. Where the prescribing health care professional has determined that an interruption in the supply of oxygen, for any reason, may have serious consequences to the user, an alternate source of oxygen should be available for immediate use.
- 25 Geriatric or any other patient unable to communicate discomfort, or hear or see the alarms while using this device, may require additional monitoring.
- 26 Oxygen supports combustion. Oxygen should not be used while smoking or in the presence of an open flame.
- 27 Do not use the device in the presence of a flammable anaesthetic mixture in combination with oxygen or air, or in the presence of nitrous oxide.
- 28 Do not use oil or grease on the concentrator or its components as these substances, when combined with oxygen, can greatly increase the potential for a fire hazard and personal injury.
- 29 If you notice any of the following, discontinue use and contact your home care provider:

- unexplained changes in the performance of this device
  - unusual or harsh sounds
  - dropped or mishandled device or the power supply
  - water spilled into the enclosure
  - broken enclosure
- 30 Use only with OX- 2A AC power supply.
  - 31 Use only with OX- 2A batteries.
  - 32 Use only approved OX- 2A accessories.
  - 33 Repairs and adjustments must be performed by OXTM authorized service personnel only. Unauthorized service could cause injury, invalidate the warranty, or result in costly damage.
  - 34 Periodically inspect electrical cords, cables, and the power supply for damage or signs of wear. Discontinue use and replace if damaged.
  - 35 To avoid electric shock, unplug the device and remove the batteries before cleaning the bag. Do not immerse the device in any fluids.
  - 36 Your home care provider is responsible for performing appropriate preventive maintenance at the intervals recommended by OXTM.
  - 37 For proper operation, your device requires unobstructed ventilation. Always make sure any openings in the case are not obstructed by items which may impede ventilation. Do not place the device in a small closed space (such as a closet). The device should not be used adjacent to or stacked with other equipment. For more information, contact your home care provider.
  - 38 Do not use an extension cord.
  - 39 Device operation above or outside of the voltage, breath rate, temperature, humidity and/or altitude values specified may decrease oxygen concentration levels.
  - 40 Never drop or insert any object into any opening.
  - 41 Be aware that the electrical cord and/or tubing could present a tripping or strangulation hazard.
  - 42 Use only power cords supplied by OXTM for this device. Use of power cords not supplied by OXTM may cause overheating or damage to the device and may result in increased emissions or decreased immunity of the equipment or system.
  - 43 Do not operate without the battery installed and working. If primary power is lost with no battery in place, the device will stop operating without warning the user. If the device must be operated without the battery, the user should be aware that there is no backup power and no warning if primary power is lost.
  - 44 Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this manual.
  - 45 Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment. See the EMC section of this manual for distances to observe between RF Generators and the E2 device to avoid interference.
  - 46 Do not use this device while sleeping unless prescribed by your clinician and use a pulse-oximeter to monitor the SpO<sub>2</sub> of patient's.

47 Do not use this device over 14h with 95% concentration oxygen and consult with doctor before oxygen therapy.

### 3.3 Cautions Overview

- 1 Keep away from heat sources (fireplaces, radiant heaters, etc.) that could cause the operating temperature at or near the device to exceed 40°C (104°F).
- 2 Don't position the equipment so that it is difficult to operate the disconnection device.
- 3 The display may be difficult to read under bright lighting conditions (sunlight, interior lights, etc.), move away from direct light for viewing the display.
- 4 Keep away from lint or other loose material that could block the intake vents.
- 5 Some countries restrict this device to be sold by or on an order of a prescribing clinician. Please ensure you comply with relevant local laws.
- 6 Non-prescribed oxygen therapy can be hazardous under certain circumstances. Use this device only when prescribed by a clinician.
- 7 Always operate the device at the setting prescribed by a clinician. Do not alter the setting unless prescribed by a clinician. Periodic reassessment of the flow settings should be done by a clinician.
- 8 Do not use this device while sleeping unless prescribed by your clinician.
- 9 It is recommended for an alternate source of oxygen to be made available in the event of power outage or mechanical failure. Consult your home oxygen provider or clinician for an appropriate backup system.
- 10 This device may not reach specified oxygen concentration purity until it has been in use for up to 2 minutes at set flowrate.
- 11 This device is designed for use by one patient at a time.
- 12 If you are unable to hear or see alarms, do not have normal tactile sensitivity, or cannot communicate discomfort, consult a clinician before using this device.
- 13 If oxygen concentration drops below the specified level, an alarm will indicate this condition. If alarm persists, stop using this device, switch to an alternate source of oxygen, and contact your home oxygen provider.
- 14 Only use approved accessories or cannula with this device. Using unapproved accessories or cannula may impair the performance of this device.
- 15 This device is not designed for use with a nebulizer. If a nebulizer is used with this device, performance may be diminished and the device may be damaged.
- 16 Always follow cannula manufacturer's instructions for proper use.
- 17 Replace the cannula on a regular basis. Check with your home oxygen provider or clinician to determine how often the cannula should be replaced.
- 18 Do not use cleaning agents other than those specified in this manual. Allow the cleaning solution to dry from the cleaned surface before use.
- 19 Always turn off this device when not in use.
- 20 Always disconnect power and turn off this device before cleaning.
- 21 Do not obstruct air intake or exhaust vents when operating this device. Blockage can cause buildup of internal heat and shut down or damage this device.

- 22 Do not place objects on top of this device.
- 23 Keep away from children and pets to prevent damage to the device and accessories and/or inadvertent setting changes.
- 24 Keep the device away from pets and pests.
- 25 Always use in a well ventilated location.
- 26 Always follow the maintenance schedule as specified
- 27 If this device indicates an abnormal condition, see Chapter 9. Troubleshooting.
- 28 Use caution when touching this device in high ambient temperatures.
- 29 Do not immerse the device or allow any liquid to enter the enclosure.

### **3.4 Overview of Important Information**

- 1 The patient is an intended operator.
- 2 Inhale through the nose for the concentrator to work most effectively. Inhaling through the mouth may result in less effective oxygen therapy.
- 3 If you have allergic reactions to nasal cannula, please contact your physician, therapist local or home care provider for assistance.
- 4 This oxygen concentrator can operate in continuous mode or pulse mode.
- 5 Please contact with manufacturer and local competent authority if have any serious incident related with the oxygen concentrator.
- 6 Manufacturer will make available on request circuit diagrams, component part lists, descriptions, calibration instructions, or other information that will assist service personnel to repair those parts of equipment that are designated by the manufacturer as repairable by service personnel.
- 7 2 hours is required for the oxygen concentrator to warm from -20°C between uses until the oxygen concentrator is ready for its intended use when the ambient temperature is 20°C, or, to cool from 50°C between uses until the oxygen concentrator is ready for its intended use when the ambient temperature is 20 °C.
- 8 The circumstances in which the user should consult a healthcare professional:
  - a) Patients with a fast breathing rate (more than 20 breaths/min) requiring a higher oxygen setting may require more oxygen than this device can produce. This device may not be appropriate in that case. Consult your clinician for alternative treatment.
  - b) It is recommended for an alternate source of oxygen to be made available in the event of power outage or mechanical failure. Consult your home oxygen provider or clinician for an appropriate backup system.
  - c) If you are unable to hear or see alarms, do not have normal tactile sensitivity, or cannot communicate discomfort, consult a clinician before using this device.
  - d) If your symptoms, such as headaches, drowsiness, confusion, fatigue or increased irritability, are not alleviated after inhaling oxygen, please consult your physician.
  - e) If you feel discomfort or are experiencing a medical emergency while undergoing oxygen therapy, seek medical assistance immediately to avoid harm.

## **4. Instructions and Training**

The Regulation (EU) 2017/745 states that the product provider must ensure that all users of this device are provided with the user manual and are fully trained in the use of the equipment.

**⚠ WARNING:** Do not use the product without proper training! Patients and care givers must be trained by an experienced person who has been authorized by the manufacturer and has appropriate training, knowledge and experience.

For further information about training contact your home oxygen provider.

## 5. Product Description

### 5.1 Device and Accessories Description

The OX- 2A Portable Oxygen Concentrator, its features, and its accessories are described in detail in this manual. Read and understand it completely before operating the device.

This manual applies to the following accessories:

- Intake filter



- AC Power Supply



- Carrying Bag



- Pull Cart



- Battery Charger



- Rechargeable LithiumIon Battery



- DC Power Cable



Power Output Plug(Plug into Concentrator)

DC cigarette lighter power plug OR Battery plug

## 6. General Instructions Before Use

A variety of accessories can enhance the portability and use of the OX- 2A Portable Oxygen Concentrator. In addition to the device, the package contains accessories to get started and a user manual. Contact your home oxygen provider for a complete list of available accessories.

Always inspect the device and its accessories for any sign of damage before use.

**⚠ IMPORTANT :** While the box or packaging may exhibit some damage, e.g., tears or dents, the device may still be in a usable condition. If the device or any accessory shows any sign of damage, contact your home oxygen provider.

Before you get started, check to make sure you have the following:

- Concentrator
- Battery
- Carry bag
- AC power supply
- DC Power Cable
- Pull Cart

### 6.1 Accessories List

Only use power supplies/adapters or accessories specified in this manual. Using accessories that are not specified may create a hazard and/or negatively affect the performance of the device.

- AC power supply - Input :100~240VAC, 50/60Hz, 2.8A, Output 13.5VDC, 14A
- DC power cable
- Carry bag
- Pull cart

**⚠ WARNING:** Do not use the device or any accessory that shows any sign of damage.

### 6.2 Battery

OX- 2A Portable Oxygen Concentrator can always be used when directly connected to a power source. However, to enhance its portability, the concentrator is equipped with a rechargeable lithium- ion external battery.

**⚠ IMPORTANT:** Optional power cords are available for various global use and travel (see Chapter 6.1.Accessories List).

Press the "check" button on the battery it will display the remaining battery power."100%"indicator light is on,indicating that the battery is fully charged.Only "Low" indicator light is on,indicating that the battery need to charge.



Fully charged



Need to charge

Fig.1

## Charging the Battery

**⚠ CAUTION:** Only charge the battery in this device or in an approved charger. (See Chapter 6.1. Accessories List.)

The charger is universal and supports a wide variety of international markets, so it can be plugged into an outlet with 100- 240V AC, 50-60 Hz. Allow one battery to charge for a minimum of four hours before use. Once completely charged, the device can run for up to 4 hours with one battery in pulse mode, at 20 breaths per minute.

**⚠ IMPORTANT:** Battery run time may vary based on breathing rate, age of battery, and environmental conditions. The concentrator is connected to a power source. The LCD display will indicate whether the device is operating on battery or external AC power.

**⚠ IMPORTANT:** Ensure power status icon (see Fig. 4) indicates power is connected. If not, check that cord is plugged in completely. (See Chapter 9. Troubleshooting for more information.)

**⚠ IMPORTANT:** After 300 charge/discharge cycles, the battery capacity will be at least 80% of its original capacity. Replace the battery when the reduced battery life is affecting your mobility.

## 6.3 Nasal Cannula

Only use a nasal cannula with the following specifications:

- 7ft (2.1 m) long
- High flow
- Crush resistant
- Large internal diameter bore
- Straight non-tapered tips
- Suitable for up to 15 liters per minute (lpm) at a max. pressure of 3.6 psi
- Suitable for adult
- Meets substance compatibility of IEC/EN 60601- 1
- Have CE marking

Recommended: Disposable Nasal Oxygen Cannula RMK01 (Produced by Ningbo Runmai Medical Technology Co.,Ltd)

**⚠ CAUTION:** Only use approved accessories with this device. Refer to the approved accessories guide for a complete list of accessories and cannula approved for use with this device. Using unapproved accessories or cannula may impair the performance of this device, including flow rate or oxygen purity. Contact your distributor for updated information and accessories or if additional, optional, or replacement accessories are needed.

## 6.4 Humidifier

Only use a humidifier with the following specifications:

- Bottle humidifier: single person multiple use
- The pressure relief valve is 3Psi and
- The maximum water Level is about 240ml
- Oxygen flow rate can be suitable for prescription from patient.
- DISS Oxygen Inlet Connector
- Outlet Port: Tapered Outlet Accepts Universal Supply Tubing End Connector

- The material shall meet requirement of ISO 18562 series.
- Have CE marking

**⚠ CAUTION:** Only use approved accessories with this device. Refer to the approved accessories guide for a complete list of accessories and cannula approved for use with this device. Using unapproved accessories or cannula may impair the performance of this device, including flow rate or oxygen purity. Contact your distributor for updated information and accessories or if additional, optional, or replacement accessories are needed.

### 6.5 Pull Cart

When using the device with a pull cart, attach and secure the concentrator with the velcro as shown in Fig. 1. The handle can be pulled out and adjusted for comfort.

**⚠ IMPORTANT:** It is recommended that patients use the pull cart to transport the device whenever possible.



Fig.2

## 7. Operating OX-2A

**⚠ IMPORTANT:** Read Chapter 3. Safety Instructions before using this device.

OX- 2A Portable Oxygen Concentrator is designed for ease of use, with all functions accessed through just a few keys on the control panel.

The device should be carried in its carry bag, placed on a cart and used when positioned upright on a table or on the floor while in the carry bag. The patient should be within the recommended cannula length during use.

### 7.1 Before Operating Your Oxygen Concentrator

1. Before operating your unit always checks to be sure the air filter (located on the back of your unit) is existed.
2. Attach the appropriate oxygen accessories to the oxygen outlet.

#### Oxygen Tubing Connection:

- A. Attach the oxygen tubing directly to the connector.(Fig.3)



Fig.3

- B. Oxygen Tubing Connection with Humidification:

**⚠ CAUTION:** Only use humidification in continuous mode.

#### **⚠ WARNING**

- Do not overfill humidifier.
- Do not reverse the oxygen input and output connections. Water from the humidifier bottle will travel through the cannula back to the patient.
- Do not use the humidifier while the device is in Pulse mod

- A. If your physician has prescribed an oxygen humidifier as part of therapy, follow these steps:

To attach the humidifier:

- Remove the cover on the humidifier bottle according to the instructions.
- Fill the humidifier with boiled tap water which has been left to cool or distilled water .
- Fill the humidifier bottle to the level specified by the humidifier bottle manufacturer mark then replace the cover according to the instructions.
- Connect the connector and the humidification bottle inlet fitting with a tube.(Fig.4)
- Attach the oxygen tubing directly to the humidifier bottle outlet fitting.(Fig.5)

Make sure all connections are secure.

Your physician has prescribed a nasal cannula, in most cases, it is already attached to the oxygen tubing. If not, follow the manufacture’s instructions for attachment.

Remove the power cord completely from the line cord strap, Make sure the power switch is in the “Off” position, and insert the plug into the wall outlet; the unit is double insulated to guard against electric shock.



Fig.4



Fig.5

**⚠ WARNING:** Improper use of the power cord and plugs can cause a burn, fire or other electric shock hazards . Do not use the unit if the power cord is damaged.

## 7.2 Connecting Nasal Cannula

**⚠ CAUTION:** Replace the cannula on a regular basis.

Check with your home oxygen provider or clinician to determine how often the cannula should be replaced.

**⚠ CAUTION:** Always follow cannula manufacturer’s instructions for proper use. Connect the tubing to the cannula port as shown in Fig.5.

1. Without humidifier: Connect the tubing to the cannula port as shown in Fig.3 if no humidifier bottle is prescribed.
2. With humidifier: Connect the tubing to the outlet of humidifier bottle as shown in Fig.5 if humidifier bottle is prescribed.
3. To connect the cannula to the patient, position the cannula tips in patient’s nostrils and pass tubing over both ears and under chin. Follow manufacturer’s instructions. Slide adapter up tubing to adjust for comfort and fit. Shown in Fig.6

Once the cannula is secured, breathe normally through the nose. OX- 2A will detect a breath and deliver the oxygen during inhalation.

**⚠ IMPORTANT:** Improper cannula placement may result in the device being unable to detect all respiratory efforts of the patient. Ensure cannula is connected securely and it has been fully inserted.





Fig.6

### 7.3 Device Overview




Item	Description	Function
1	Control Panel	Control switches and LCD display
2	Oxygen outlet	Oxygen is dispersed through this port
3	Humidifier Holder	To hold the humidifier
4	Air Intake Vent	Air inlet for enclosure ventilation fan
5	Carrying Handle	Hand grip area for transporting the device
6	Filter cover	For filter replacement
7	Air Exhaust Vent	Air outlet for enclosure ventilation fan
8	Power Input Connector	Connection point for OX-2A external power supplies: AC adapter, vehicle DC, external Battery
9	Power Switch	I =ON , O=OFF



## 7.4 Control Panel




Item	Description	Function
1	Power	Turns the device On and Off
2	Mode	Selects one of the two possible operating states of the device
3	Display Screen	Shows information about the operating status of the device
4	Green light	Illuminates when your concentrator is operating and the oxygen concentration $\geq 82\%$
5	Yellow light	The oxygen concentration $< 82\%$ and indicates caution or attention required
6	Plus/Minus	Increases the displayed setting in cycle
7	Bell, cancel temporary	Cancel the bell temporarily and/or reuse it


## 7.5 Run Time Screen Description




Item	Description	Function
1	External Power Status	Indicate external power supply
2	Operating Mode	 :Continuous Mode,  :Pulse Mode
3	Alarm Status	Indicate the type of alarm: NF:No Flow Alarm NB:No Breath Alarm LO:Low Oxygen Concentration Alarm HB:High Breath Rate Alarm LB:Low Battery Alarm
4	Setting	Indicate the flow rate setting
5	Battery Status	Indicate battery supply and its quantity
6	Total Running times	Indicate the usage time
7	Audible Alarm Status	Indicate that it is in the ringing alarm state
8	Breath Indicator	Indicate breath detected in pulse mode


## 7.6 Turning On

- To turn the device on, Press the power switch to the “On” position, Then press the power key .
- The concentrator will chirp and the green, yellow LEDs will flash once,  
Yellow LED - indicates caution or attention required  
Green LED - indicates device is on. The green LED will then stay lit.

 **IMPORTANT:** No adjustments can be made until the startup sequence is completed.


## 7.7 Adjusting Setting

 **IMPORTANT:** After powering on OX- 2A, the startup sequence will take approximately 35 seconds. Specified oxygen level will be reached within 2 minutes of use.

 **IMPORTANT:**The device is for prescription use, the quantity of concentration of Oxygen(flow rate and use time) to delivered shall be set according to prescription from your physician.


- The device starts working in the previous setting.

### Mode switch:



- Use the mode button  to alternate between Pulse mode and continuous mode as Fig.7 showed.
- In pulse mode, the device will deliver a pulse of oxygen at the beginning of each of your inhalation.
- In continuous mode, the device will provide a continuous flow of oxygen, but will consume more power and have a shorter battery life.

### Setting the mode can be done as follows:



**Note:** When it is turned on, the device will automatically start at the flow rate setting used when the device was last turned off. As a precaution, each time you start the device, verify the flow setting.




- To set the device's flow setting, press the  button.

### Pluse Mode:

- The device can be set to 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11 and 12 in Pulse Mode. Press the  button to increase the flow rate from 1 to 12 in 1 increment. When the setting is increased to be 12, press the  button once again, the setting will return to 1.

**Continuous Mode:**

- The device can be set to 0.5, 1, 1.5, and 2 in Continuous Mode. Press the  button to increase the flow rate from 0.5 to 2 in 0.5 increment. When the setting is increased to be 2, press the  button once again, the setting will return to 0.5.

-  **IMPORTANT:** If an air leak is suspected, leaks can be detected with a solution of soap and water applied to the cannula-concentrator connection point and looking for bubbles.
-  **WARNING:** It is very important to set your device to your prescribed level of oxygen flow. Do not increase or decrease your flow rate from your prescribed level until you first consult with your physician.
-  **IMPORTANT:** Flow can be verified by putting the oxygen concentrator in continuous mode and placing the end of the nasal cannula under the surface of a half full cup of water and looking for bubbles. The current setting and power source (external power or battery; battery icon also shows approximate






Fig.7


**7.8 Responding to Alarms**

-  **CAUTION:** If you are unable to hear or see alarms, do not have normal tactile sensitivity, or cannot communicate discomfort, consult your clinician before using this device.
-  **IMPORTANT:** The alarm system is tested during the startup sequence. You should see all alarm lights briefly turn on and the audible alarm indicator chirp. If alarms are suspected of miss-operating, contact your distributor for verification that alarms are working correctly.

## 7.9 Turning Off

 **CAUTION:** Always turn off this device when not in use.




To turn the OX- 2A Portable Oxygen Concentrator off, press and hold the power key. The device will chirp and the screen will display a shutdown message for approximately five seconds, then go into low-power mode.




 **IMPORTANT:** Do not disconnect the AC power supply at the same time while the unit is running. Always use the power key to turn the device off. Wait until the device has completely shut down before disconnecting from power.

## 8. Alarm

OX- 2A oxygen concentrator has 2 types of alarms according to different potential risks during the alarm active: Medium Priority and Low Priority.

When alarm occurs, it will remind the user by the alarm indicator light, alarm sound and alarm messages, different types of alarms have different alarm information as follows:

Alarm Code	Visual, Auditory	Description	Priority Alert	Alarm Condition	Alarm delay	Alarm Signal	Remedy
LO	 <p>3 short beep every 11.2 seconds + yellow light flashes</p>	<p><b>Low Oxygen Concentration Alarm</b></p> <p>This alarm occurs when the device is delivering a lower concentration of oxygen than specified(82%). Operator shall make instant on-time response.</p>	Medium	Technical Alarm	<5s	Visual + Auditory	<ol style="list-style-type: none"> <li>1.Change to another source of oxygen and contact your equipment provider.</li> <li>2. Filter restricted: Clean or replace air inlet filter. Place your device so there is adequate air flow.</li> <li>3.Restriction in humidifier or tubing: Repair or replace humidifier or tubing.</li> </ol>
NB	 <p>2 short beep every 16 seconds + yellow light on</p>	<p><b>No Breath Alarm</b></p> <p>This alarm occurs when a breath is not detected for a period of 2 minutes or more. This alarm becomes silent as soon as a breath is detected. If no breath is detected after approximately 30 minutes, the device shuts down to conserve power.Operator shall pay attention.</p>	Low	Technical Alarm	<5s	Visual + Auditory	Check the connection from the cannula to the device. Make sure the nasal cannula is properly positioned on your face and that you are breathing through your nose. Make sure the cannula tubing is not kinked or obstructed.
HB	 <p>2 short beep every 16 seconds + yellow light on</p>	<p><b>High Breath Rate Alarm</b></p> <p>This alarm indicates that the user's breath rate is exceeding the capacity of the device. The device is still working properly and is still providing oxygen at the rate specified for the setting.Operator shall pay attention.</p>	Low	Technical Alarm	<5s	Visual + Auditory	The indicator resets itself when the breath rate is over 40 time. If this indicator regularly occurs, contact your home care provider.

Alarm Code	Visual, Auditory	Description	Priority Alert	Alarm Condition	Alarm delay	Alarm Signal	Remedy
LB	 <p>2 short beep every 16 seconds + yellow light on</p>	<p><b>Low Battery Alarm</b></p> <p>This alarm occurs when approximately 30 minutes of battery life remains. Remaining battery life is dependent on your device settings and your activity level. The empty battery symbol flashes on-screen. Operator shall pay attention.</p>	Low	Technical Alarm	<5s	Visual + Auditory	Replace the battery or connect to a power source.
NF	 <p>2 short beep every 16 seconds + yellow light on</p>	<p><b>No Flow Alarm</b></p> <p>This alarm occurs when the device detects there is no oxygen flowing in the patient cannula. Operator shall pay attention.</p>	Low	Technical Alarm	<5s	Visual + Auditory	<ol style="list-style-type: none"> <li>1. Check the cannula for kinks or other obstructions that are stopping oxygen flow through it.</li> <li>2. Check whether Filter Blocked or not, if so, clean or replace air inlet filter.</li> </ol>
ON	 <p>2 short beep every 16 seconds + yellow light on</p>	<p><b>Oxygen flow outside normal limits alarm</b></p> <p>This alarm occurs when the device detects there is oxygen flowing outside normal limits in the patient cannula. Operator shall pay attention.</p>	Low	Technical Alarm	<5s	Visual + Auditory	<ol style="list-style-type: none"> <li>1. Check the cannula for kinks or other obstructions that are stopping oxygen flow through it.</li> <li>2. Check whether Filter Blocked or not, if so, clean or replace air inlet filter.</li> <li>3. If not issues mentioned above, change to another source of oxygen and contact your equipment provider.</li> </ol>
	<p>(Power Connection Symbol Disappears)</p> <p>2 short beep every 16 seconds + yellow light on</p>	<p>The Power Connection symbol disappears when the device is disconnected from the power source. Operator shall pay attention.</p>	Low	Technical Alarm	<5s	Visual + Auditory	Check whether the power supply is connected

Sound pressure level range of the auditory is more than 65dB.

**⚠ IMPORTANT:** After running 5 minutes, disconnect the device from the power supply by unplugging the power cord, an audible alert sounds like “beep” and the yellow light on, which means the alarming system of the device works normally.

**⚠ IMPORTANT:**All alarms have been set and cannot be changed.

**⚠ CAUTION:** The operator needs to operate the oxygen concentrator in front of the control panel of the oxygen concentrator.

**⚠ CAUTION:** When the power failure, power failure alarm will occur and other alarms will disappear. If power is lost for less than or equal to 30s, power failure alarm will disappear and the alarm prior to the power lost will occur again.

## 9. Troubleshooting

**IMPORTANT:** The table below lists common problems and actions you can take. If you are unable to resolve a problem, please contact your home care provider.

Problem	Possible Cause	What You Should Do
Device won't turn on	Battery is depleted.	Use the AC or DC power cords to operate the device. If this does not resolve the problem, contact your home care provider.
	AC power connection is not correct	Re-install it correctly.
	No Battery connected	Connect battery or external power
Device will not trigger a pulse of oxygen	OX- 2A is not turned on.	Press the power button.
	Tubing/cannula longer than 7 feet (2.1m)	Attach 7 feet (2.1m) tubing/cannula
	Cannula tubing is kinked or twisted.	Make sure the tubing is connected properly to the oxygen outlet port and that it is free of any obstruction.
	Humidifier attached	Remove humidifier
	No inspiration detected	Contact your dealer/provider
	Unit in Continuous Mode	Switch to Pulse Mode
	Device malfunction.	Contact your home care provider.
Oxygen not at full concentration	Device is warming up.	Wait 5 minutes for the device to deliver oxygen at the prescribed concentration. If the condition persists, contact your home care provider.
	Filter blocked	Clean or replace air inlet filter. Place your device so there is adequate air flow.
	Dirty, or Blocked, or defective humidifier bottle.	Detach the humidifier from the oxygen outlet. If proper flow is obtained, clean or replace humidifier.
	Blocked or defective cannula, or oxygen tubing.	Detach cannula. If proper flow is restored, clean or replace if necessary. Disconnect the oxygen tubing at the oxygen outlet. If proper flow is restored, check oxygen tubing for obstructions or kinks or defective. Replace if necessary.
	compressor failure	Contact your dealer/provider
	Sieve beds may be faulty or use up	Contact your dealer/provider
	Gas leakage	Check connection between the device outlet and tubing, connection between the tubing and humidifier bottle, have any gas leakage. Replace if necessary.
	Device malfunction.	Contact your home care provider.
Low oxygen flow	Filter restricted:	Clean or replace air inlet filter. Place your device so there is adequate air flow.
	Dirty, or Blocked, or defective humidifier bottle.	Detach the humidifier from the oxygen outlet. If proper flow is obtained, clean or replace

Problem	Possible Cause	What You Should Do
		humidifier.
	Blocked or defective cannula, or oxygen tubing.	Detach cannula. If proper flow is restored, clean or replace if necessary. Disconnect the oxygen tubing at the oxygen outlet. If proper flow is restored, check oxygen tubing for obstructions or kinks or defective. Replace if necessary.
	Gas leakage	Check connection between the device outlet and tubing, connection between the tubing and humidifier bottle, have any gas leakage. Replace if necessary.
	Other	Contact your dealer/provider
Alarm Occurs	The device needs your attention.	See the Alarm Indicators and Screen Symbols section for information about specific alarms and what you should do.

## 10. Maintenance and Cleaning

### 10.1 Routine Maintenance

**⚠ WARNING:** Do not use lubricants on this device or any of its accessories.

**⚠ CAUTION:** Replace the cannula on a regular basis. Check with your distributor or clinician to determine how often the cannula should be replaced.

Device will indicate with an alarm when a filter or component needs to be leaned or changed. (Also, see Chapter 9. Troubleshooting.)

**⚠ IMPORTANT:** The cannula and patient filter can be contaminated from the patient, care in handling these components should be taken.

#### Air Filter and Oxygen Outlet Connector

The connector should not be cleaned at least once a week. To clean, follow these steps:

1. Remove the oxygen outlet connector (if used ).
2. Wash in a solution of warm water and dishwashing detergent.
3. Rinse thoroughly with warm tap water and towel dry .The filter should be completely dry before reinstalling.

Air filter need to be replaced after device operates each 1500 hours. Air filter is disposable after used. Not replace the air filter regularly will affect the performance of concentration.

**⚠ CAUTION:** To prevent product damage, do not attempt to operate the unit without the air filter or while the filter is still damp.

### 10.2 Cleaning

**⚠ WARNING:** Do not submerge this device in liquid. Do not expose to water or precipitation. Do not expose to dusty conditions.

**⚠ CAUTION:** Do not use cleaning agents other than those specified in this manual. Allow the cleaning solution to dry from the cleaned surface before use.

**⚠ CAUTION:** Always disconnect power and turn off this device before cleaning.

Clean the exterior with a soft cloth slightly dampened with soapy water or with anti-bacterial wipes (Isopropyl alcohol 70% solution).

**▲ IMPORTANT:** The device should receive an external cleaning weekly, accessories should be cleaned as needed. And the device can be cleaned by bacterial wipes over 156 cleaning cycles.

Nasal cannula: Refer to the original manufacturer's instructions for cleaning the nasal cannula.

### **10.3 Service Life**

The expected service life of the device is 3 years or 10,000 hours, except for the sieve beds. The service life of the sieve beds will depend on the operating conditions. Replace them as needed, indicated by the low oxygen concentration alarm. If intake and exhaust vents are not blocked and the low oxygen concentration alarm persists, contact your distributor for instructions on replacing the sieve beds.

Air Filter replacement time: 1500 hours.

## **11. Device Repair and Disposal**

### **11.1 Repair**

Do not attempt to repair the device. Contact your home oxygen provider or distributor for assistance (see Chapter 9. Troubleshooting).

### **11.2 Disposal**

This device contains electrical and/or electronic components that must be recycled per EU Directive 2012/19/EU-Waste Electrical and Electronic Equipment (WEEE).

- Contact your distributor regarding disposal of the device.
- Battery is rechargeable and can be recycled, always properly dispose of battery according to local regulations or contact your distributor.

## **12. Warranty**

The device warranty is limited to three (3) years from date of manufacture or 10,000 hours of total use. All accessories including the batteries are limited to two (2) year warranty.

The standard warranty is only valid for products handled as stated in the user manual and in accordance with general industry good practice and standards.

## **13. Disclaimer**

### **13.1 Disclaimer**

The information in this document has been carefully checked and is believed to be reliable.

The manufacturer reserves the right to modify any of these products to enhance readability, functionality or design. The manufacturer does not assume any liability arising from the application or use of any product or circuit described. It does not include any license under its patent rights or the rights of others.

### **13.2 This Document**

The information in this document is subject to change without any notice. This document contains proprietary information for copyright protection. This document shall not be reproduced in whole or in any form without the prior written consent of the manufacturer (except for a brief excerpt of the reviews and scientific papers). Be sure to read and understand all the manuals provided by the product.

For help

Please contact your domestic oxygen supplier or distributor if you have any concerns about the information in this manual or the safety operation of the device.

## 14. Specifications

Operating Conditions	Operational temperature: 50°F to 104°F (10°C to 40°C) Relative humidity: 15% to 95% Altitude: up to 7,500 ft. (2500m.) Atmospheric pressure 50-106 kpa (7.3-15.4 psi)
Storage and Transport Conditions	-4°F to 122°F (- 20°C to 50°C) - unit only Relative humidity: up to 93%, non-condensing Atmospheric pressure 50-106 kpa (7.3-15.4 psi)
Power supply	External power supply: 230V~, 50Hz, 0.7A Battery supply:DC14.8V,10A
Maximum rated voltage	240V
Backpressure	7Kpa if humidifier bottle attached 0 Kpa if no humidifier bottle attached
Oxygen Concentration	90% ± 3% at all settings
Flow Settings and Pulse Volumes	Pulse Mode 1 = 6 ml; 2 = 12 ml; 3 = 18 ml; 4 = 24 ml; 5 = 30 ml; 6 = 36 ml; 7 = 42 ml; 8 = 48 ml; 9 = 54 ml; 10 = 60 ml; 11 = 66 ml; and 12= 72 ml; +/- 15% or 4 ml, whichever is greater (Average of 20 consecutive breaths) Inspiratory trigger pressure: adjustable between -0.125cmH2O to -0.4cmH2O up to a max of 2000 ml/min +/- 200 ml Continuous Mode 0.5 = 500 ml/min; 1 = 1000 ml/min; 1.5 = 1500 ml/min; 2 = 2000 ml/min; +/- 10% or 100 ml/min, whichever is greater (3 minute running average) Note: Max recommended flow is 2 LPM (at nominal outlet pressures of 0 and 7 kPa).
User Interface	Power switch, Push buttons, back- light liquid crystal display (LCD)
Dimensions	200 (L) X 230(W) X 325(H)mm
Weight	5.4Kg
Sound Level	46 dBA max. at Nominal Continuous Mode Setting of 2(1M from front)
Audio Alarm	50 decibels (nominal) at 39 in. (1M from front)
Outlet Pressure	8.5 psi +/- 0.5psi
Equipment Class and Type	Degree of Protection Against Electric Shock: Class II Degree of Protection Against Electric Shock: Type BF Applied Part (Nasal Cannula) Degree of Protection Against ingress of water and particulate matter: IP22 Mode of Operation: Continuous/Pulse
Sieve Bed	Dimensions:148 (L) X 71(W) X 214(H)mm Weight:1kg ± 1%

	Outlet Pressure:8.5 psi $\pm$ 0.5psi Outlet Size: $\phi$ 7 Molecular Sieves: Specification/type:NITROXY-REVOLUTION Manufacturer:CECA SA
Total work hours	10,000 hours

**Pulse mode bolus size (ml/breath) versus setting and breath rate**

Breath per minute	Setting											
	1	2	3	4	5	6	7	8	9	10	11	12
15	6	12	18	24	30	36	42	48	54	60	66	72
20	6	12	18	24	30	36	42	48	54	60	66	72
25	6	12	18	24	30	36	42	48	54	60	66	72
30	6	12	18	24	30	36	42	48	54	60	66	72
35	6	12	18	24	30	36	42	48	54	60		
40	6	12	18	24	30	36	42	48	54			

All values +/- 15% over all operating conditions

**Continuous Mode Flow (L/min) versus setting**

Setting	Flow rate
0.5	0.5
1.0	1.0
1.5	1.5
2.0	2.0

All values +/- 0.2 liters over all operating conditions

**Battery duration in continuous mode**

Setting	Minutes
0.5	240
1.0	165
1.5	120
2.0	75

*Battery duration times will degrade with battery age,use over time,and operational condition of the POC*

**Battery duration in pulse mode(20BPM)**

Setting	Bolus size (ml/breath)	Minutes
1	6	300
2	12	235
3	18	235
4	24	200
5	30	200
6	36	165

7	42	165
8	48	145
9	54	135
10	60	120
11	66	115
12	72	110

**Battery duration times will degrade with battery age, use over time, and operational condition of the POC**

**⚠ CAUTION:** The battery needs to be replaced when the low battery alarm is occurred after fully charged battery is connected.

## 15. Traveling With Your OX-2A

With your OX- 2A and proper advance planning, you can enjoy traveling within your community and beyond. Before you leave, make sure you pack the following:

- Fully charged battery (and extra batteries for a long trip).
- AC power supply and connector cord
- DC power supply
- Carrying bag
- Mobile cart

Also, be sure to take the telephone numbers of your home care provider and physician in case of an emergency.

### 15.1 By Motor Vehicle


















Use the OX- 2A DC power supply to plug in your device using the vehicle’s cigarette lighter or DC power input. You can also use the OX- 2A device while it is running from a DC power source.




**⚠ CAUTION:** Make sure the vehicle has been started before plugging in your DC power supply. If you operate the device using the DC power cord while the vehicle engine is turned off, you could inadvertently drain the vehicle’s battery.

### 15.2 By Bus or Train

Most bus and train lines allow passengers to use portable oxygen concentrators, but you may need to notify them in advance. When you make your travel arrangements, contact your carrier well before your departure for permission to bring your system and use it on board.

## 16. Symbols

	Follow instructions for use
	Type BF applied part
	Class II equipment
	Bell, cancel temporary
	Stand-by button
	Machinery Directive
	Manufacture's serial number
	Authorized representative in the European community
	Date of manufacture
	Manufacture
	This device contains electrical and/or electronic components that must be recycled per EU Directive 2012/19/EU-Waste Electrical and Electronic Equipment (WEEE).
	The equipment bears CE mark CE 2460 indicating its conformity with the provision of Regulation(EU) 2017/745 concerning medical devices.
	FRAGILE Contents of the transport package are fragile therefore it shall be handled with care.
	THIS WAY UP Indicates correct upright position of the transport package.
	KEEP AWAY FROM RAIN Transport package shall be kept away from rain.
	Non-ionizing electromagnetic radiation
	Enclosure protection classification "2" means protection against solid foreign objects of $\phi 12.5$ mm and greater "2" means protection against falling water dripping (15° tilted)

	No sitting
	No smoking
	No open flame: Fire, open ignition source and smoking prohibited

## 17. EMC Information

- ⚠ Warning:** Don't use the oxygen concentrator (model: OX-2A) near high-frequency surgical equipment, near magnetic resonance imaging equipment, or where the intensity of electromagnetic disturbances may be high.
- ⚠ Warning:** Use of oxygen concentrator (model: OX-2A) near other equipment should be avoided because it could result in improper operation. If such use is necessary, oxygen concentrator (model: OX-2A) and the other equipment should be observed to verify that they are operating normally.
- ⚠ Warning:** Use of accessories, transducers, and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of oxygen concentrator (model: OX-2A) and result in improper operation.
- ⚠ Warning:** Portable and mobile RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 12 in (30 cm) to any part of the oxygen concentrator (model: OX-2A), including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Guidance and Manufacturer's Declaration – Electromagnetic Emission		
OX- 2A is intended for use in the electromagnetic environment specified below. The user of OX- 2A should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	The Device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The Device is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3- 2	Not applicable	
Voltage Fluctuations/Flicker Emissions IEC 61000-3- 3	Not applicable	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity		
OX- 2A is intended for use in the electromagnetic environment specified below. The customer or user of OX- 2A should assure that it is used in such an environment.		
Anti-interference detection	IEC 60601 Test Level	Compliance Level
Electrostatic Discharge (ESD) IEC 61000-4- 2	Contact: +8 KV Air: +2,+4,+8,+15 KV	Same as the left
Electrical Fast Transient/Burst IEC 61000-4-4	The input a.c. power ports: ±2 KV The input d.c. power ports: ±2 KV Signal input/output ports: ±1 KV	Not applicable
Surge IEC 61000-4- 5	Input power ports: +0.5, +1.0 KV Signal input/output:+2.0 KV	Not applicable
Voltage dips IEC 61000-4-11	0.5 cycles for > 95% (sync angle (degrees):0, 45, 90, 135, 180,225, 270, 315)  1 cycles for >95% UT (sync angle (degrees):0)	Not applicable

	25 (50Hz)/30 (60Hz) cycles for 30% U <sub>T</sub> (sync angle (degrees):0)	
Voltage interruption IEC 61000-4-11	250 (50Hz)/300 (60Hz) cycles for >95% UT (sync angle (degrees):0)	
Power Frequency (50/60 Hz) Magnetic Field IEC 61000-4-8	30A/m	Same as the left
Note: U <sub>T</sub> is the A.C. mains voltage prior to application of the test level.		

Guidance and Manufacturer's Declaration - Electromagnetic Immunity								
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band <sup>a)</sup> (MHz)	Service <sup>a)</sup>	Modulation <sup>b)</sup>	Maximum Power(W)	Distance (m)	IMMUNITY TEST LEVEL (V/M)	
	385	380-390	TETRA 400	Pulse modulation <sup>b)</sup> 18 Hz	1.8	0.3	27	
	450	430-470	GMRS 460 FRS 460	FM ±5 kHz deviation 1 kHz sine	2	0.3	28	
	710	704-787	LTE Band 13,17	Pulse modulation <sup>b)</sup> 217 Hz	0.2	0.3	9	
	745							
	780							
	810	800-960	GSM800/900, TETRA800, CDMA850, LTE Band 5	Pulse modulation <sup>b)</sup> 18 Hz	2	0.3	28	
	870							
	930							
	1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3,4,25; UMTS	Pulse modulation <sup>b)</sup> 217 Hz	2	0.3	28	
	1845							
	1970							
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID2450, LTE Band 7	Pulse modulation <sup>b)</sup> 217 Hz	2	0.3	28	
	5240	5100-5800	WLAN 802.11 a/n	Pulse modulation <sup>b)</sup> 217 Hz	0.2	0.3	9	
5500								
5785								

NOTE:  
1. If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.  
a) For some services, only the uplink frequencies are included.  
b) The carrier shall be modulated using a 50 % duty cycle square wave signal.  
c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

## 18. Declaration of Conformity



**Oxytek Medical Technology Co., Ltd.**

**ADDRESS:** 10-2 , Guang Long Industrial Park, South 1st Rd, Chencun town, Shunde district, Foshan city, Guang Dong Province , China.

**Product Designation:** Oxygen Concentrator


**Catalog Number:** OX-2A

**Tel:** +86-757-2331-1740

**Fax :**+86-757-2331- 1745

**Web:** www.oxtm-o2.cn

**Email :** jenny@oxtm-o2.cn



**WellKang Ltd (www.CE-marking.eu)**

**Enterprise Hub, NW Business Complex,**

**1 Beraghmore Road,Derry, BT48 8SE,**

**Northern Ireland**

**Web:** www.CE-Marking.eu & www.Wellkang.Ltd.uk

**Document No.:** OX-2A-002-017

**Version:** A5

**Date:** 2023-09-25

**Software version:** OX2A-1.00



## 19. Reporting System of Serious Events

Please contact with manufacturer and local competent authority as below table once any serious incident related with the oxygen concentrator occurs.

Country	Contact information	Website
Belgium	<p><b>MDD AIMDD - IVDMD</b> Head of Vigilance Division: Th. Roisin</p> <p><b>MDD AIMDD Vigilance: C. Driesmans</b> tel: +32 2 528 4418</p> <p><b>IVDMD Vigilance: J. Poels</b> tel: +32 2 528 4449 E-mail: vigilance.meddev@fagg-afmps.be</p> <p><b>FAMHP- Federal Agency for Medicines and Health Products</b> Place Victor Horta 40, box 40, B - 1060, Brussels, fax:+32 2 528 4120</p>	<a href="https://www.afmps.be/fr">https://www.afmps.be/fr</a>
Bulgaria	<p>Executive Director of BDA: <b>Bogdan Kirilov, M.Pharm.</b> Head of Division 'Medical devices': <b>Todor Darakchiev</b> <b>Bulgarian Drug Agency</b> 8 Damyan Gruev Str., BG - 1303 Sofia, tel:+359 2 890 34 83, fax:+359 2 890 34 34 E-mail: todor.darakchiev@bda.bg, bda@bda.bg - Web site</p>	<a href="https://www.bda.bg/bg/">https://www.bda.bg/bg/</a>
Ceska Republika / Czech Republic	<p><b>Ivana Justová</b> <b>State Institute for Drug Control</b> Šrobárova 48, 100 41 Prague 10, Czech Republic, tel: +420 272 185 794, fax: +420 272 185 764 E-mail: urgent@sukl.cz , ivana.justova@sukl.cz</p>	/
Hrvatska / Croatia	<p>Krunoslav Kranjcec, Agency for Medicinal products and medical devices Ksaverska cesta 4, 10 000 Zagreb, tel: +385 1 4884 327, fax: +381 1 4884 110 E-mail: medpro@halmed.hr, Krunoslav.kranjcec @halmed.hr</p>	<a href="https://www.halmed.hr/?ln=en">https://www.halmed.hr/?ln=en</a>
Danmark / Denmark	<p><b>Danish Medicines Agency</b> Axel Heides Gade 1, DK - 2300 - Kobenhavn, tel:+45 44 88 9595, fax:+45 44 88 9599 E-mail: med-udstyr@dkma.dk, Web sites: www.medicinskudstyr.dk</p>	<a href="https://laegemiddelstyrelsen.dk/en/devices/">https://laegemiddelstyrelsen.dk/en/devices/</a>
Deutschland / Germany	<p><b>AIMDD, MDD - Dr. Ekkehard Stöblein</b> - tel:+49 228 207 5384</p> <p><b>IVDMD - Prof. Dr. Rüdiger Siekmeier</b> - tel:+49 228 207 5360</p> <p><b>Federal Institute for Drugs and Medical Devices</b> Kurt Georg Kiesinger Allee 3, D - 53175 Bonn, fax:+49 228 207 5300 E-mail: medizinprodukte@bfarm.de</p>	<p><a href="https://www.bfarm.de/DE/Home/_node.html">https://www.bfarm.de/DE/Home/_node.html</a></p> <p><a href="https://www.pei.de/DE/home/home-node.html">https://www.pei.de/DE/home/home-node.html</a></p>

Country	Contact information	Website
	<b>IVDMD</b> <b>Dr. Markus Funk</b> - tel:+49 6103 77 3115 <b>Jochen Halbauer</b> - tel:+49 6103 +77 3114 <b>Paul Ehrlich Institute</b> , Section Pharmacovigilance 2 Paul-Ehrlich-Strasse 51-59, D - 63225 Langen, fax:+49 6103 77 1268 E-mail: pharmacovigilance2@pei.de	
Eesti / Estonia	<b>Sofia Ratusnaja</b> - tel:+372 744 7425 <b>Health Board, Medical Devices Department</b> 1a Põllu st., EE - Tartu 50303 E-mail: mso@terviseamet.ee -	<a href="https://www.terviseamet.ee/en/medical-devices">https://www.terviseamet.ee/en/medical-devices</a>
Ireland / Eire	<b>Health Products Regulatory Authority</b> Kevin O'Malley House, Earlsfort Centre, Earlsfort Terrace, IE - Dublin 2 E-mail: devicesafety@hpra.ie	<a href="https://www.hpra.ie/">https://www.hpra.ie/</a>
Ellada / Greece	<b>Eleni Papaioannou, MD</b> - tel:+30 213 20 40542, fax:+30 210 65 49585 E-mail: vigilancematerial@eof.gr <b>National</b> <b>Organization for Medicines</b> 284 Mesogion Ave, GR- 15562 Holargos, Athens	/
España / Spain	<b>Carmen Abad</b> <b>Carmen Valls</b> - tel: +34 91 822 5255, fax: +34 91 822 5289 <b>Agencia Española de Medicamentos y</b> <b>Productos Sanitarios</b> C/ Campezo 1, Edificio 8, ES - 28022 Madrid E-mail: psvigilancia@aemps.es	<a href="https://www.aemps.gob.es/">https://www.aemps.gob.es/</a>
France	<b>Emilie Alliez</b> - tel:+33 1 55 87 33 46, fax:+33 1 55 87 37 02 <b>Agence nationale de sécurité du</b> <b>médicament et des produits de santé</b> <b>(ANSM)</b> 143-147 boulevard Anatole France, FR - 93285 Saint Denis Cedex E-mails: Exclusively for correspondence between authorities:medicaldevicesvigilance@ansm.s ante.fr Other purposes: materiovigilance@ansm.sante.fr	<a href="http://www.afssaps.fr/">http://www.afssaps.fr/</a>
Italia / Italy	<b>Vigilance on Medical Devices</b> <b>Head of Unit 5 - Dr.ssa Lucia Lispi</b> - tel:+39 06 5994 2055 E-mail: dgfdm@postacert.sanita.it, vigilance@sanita.it, l.lispi@sanita.it  <b>MDD AIMDD Vigilance - Dr.ssa Antonella</b> <b>Campanale - Dr.ssa Daniela Minella</b> tel: +39 06 5994 3038, +39 06 5994 3069 E-mail: dgfdm@postacert.sanita.it, vigilance@sanita.it, a.campanale@sanita.it; d.minella@sanita.it	/

Country	Contact information	Website
	<p><b>Head of Unit 4 - Dr.ssa Antonella Colliardo - tel:+39 06 59943968.IVDMD Vigilance - Dr.ssa Maria Gabriella Cividino - Dr.ssa Maria Elena Russo</b>  tel: +39 06 59943785, +39 06 59942516  E-mail: dgfdm@postacert.sanita.it, mg.cividino@sanita.it; me.russo@sanita.it; a.colliardo@sanita.it  Ministry of Health, Directorate General of Medical Devices and Pharmaceutical Services Via Giorgio Ribotta 5, IT - 00144 Roma</p>	
Kypros / Kibris / Cyprus	<p><b>Ioannis Argyropoulos – tel: +357 22 605785 Cyprus Medical Devices Competent Authority</b>  Prodromou 1 &amp; Chilonos 17 Corner, CY - 1449 Nicosia, fax:+357 22 468427  E-mail: cymda@mphs.moh.gov.cy</p>	/
Latvija / Latvia	<p>Medical Device Evaluation Department - tel: +371 67 078 466, tel: +371 67078466  State Agency of Medicines, 15 Jersikas street, LV - 1003 Riga  E-mail: info@zva.gov.lv</p>	/
Lietuva / Lithuania	<p>Director: <b>Nora Ribokiene</b> - tel:+370 5 261 51 77, fax:+370 5 212 73 10  <b>The State Health Care Accreditation Agency, under the Ministry of Health of the Republic of Lithuania</b>  Jeruzales str. 21, LT-08420 Vilnius  E-mail: vaspvt@vaspvt.gov.lt</p>	<a href="https://vaspvt.gov.lt/">https://vaspvt.gov.lt/</a>
Luxembourg	<p>Ministère de la Santé, Direction de la Santé - tel : +352 247 85612  Villa Louvigny - allée Marconi, L - 2120 Luxembourg  E-mail: meddevices.vigilance@ms.etat.lu</p>	<a href="https://sante.public.lu/fr/politique-sante/ministere-sante/index.html">https://sante.public.lu/fr/politique-sante/ministere-sante/index.html</a>
Malta	<p><b>Malta Medicines Authority - Medical Devices Unit Medicines Authority</b>  Sir Temi Żammit Buildings, Malta Life Sciences Park  San Ġwann SĠN 3000, Malta  Tel : +356 2343 9000  E-mail : devices.medicinesauthority@gov.mt</p>	<a href="https://medicinesauthority.gov.mt/medicaldevices">https://medicinesauthority.gov.mt/medicaldevices</a>
Magyarország / Hungary	<p><b>Kornel Szerdi dr.- tel:+36 1 886 9329, fax:+36 1 269 1255 Health Registration and Training Centre, Department of Medical Devices</b>  1051, Budapest, Zrínyi street 3, Hungary  E-mail: amd.vig@ogyei.gov.hu</p>	/
Nederland / Netherlands	<p><b>Sietske Eerens, Esther Klinckenberg</b> - tel:+31 88 120 5000, fax:+31 88 120 5001  <b>Dutch Health and Youth Care Inspectorate, IGJ Information Office (Meldpunt)</b>  Visitors address: Stadsplateau 1   3521 AZ   Utrecht, postal address: Postbus 2518   6401 DA  </p>	<a href="https://www.igj.nl/">https://www.igj.nl/</a>

Country	Contact information	Website
	Heerlen E-mail: meldpunt@igj.nl	
Österreich / Austria	<b>Federal Office for Safety in Healthcare - (BASG) Bundesamt für Sicherheit im Gesundheitswesen</b> <b>Institute for Surveillance, Department Medical Devices Surveillance</b> Traisengasse 5, A-1200 Vienna, fax: +43 50555 36409 E-mail: medizinprodukte@basg.gv.at	<a href="https://www.urpl.gov.pl/pl">https://www.urpl.gov.pl/pl</a>
Polska/ Poland	<b>Competent Authority</b> <b>Andrzej Karczewicz</b> - tel:+48 22 492 11 90, <b>Beata Koziozemska</b> - tel: +48 22 492 11 68 <b>Office for Registration of Medicinal Products, Medical Devices and Biocidal Products</b> Al. Jerozolimskie 181C, 02-222 Warsaw, fax:+48 22 492 11 99 E-mail: incydenty@urpl.gov.pl	<a href="https://www.urpl.gov.pl/pl">https://www.urpl.gov.pl/pl</a>
Portugal	<b>Raquel Alves</b> - tel: + 351 21 798 7297, tel:+351 21 798 7145, fax:+351 211 117 559 <b>Infarmed - National Authority of Medicines and Health Products, IP Unidade de Vigilância de Produtos de Saúde da Direção de Produtos de Saúde</b> Parque da Saúde de Lisboa, Av. do Brasil, nº 53, PT - 1749-004 Lisboa E-mail: dvps@infarmed.pt	<a href="https://www.infarmed.pt/web/infarmed/entidades/dispositivos-medicos/vigilancia-de-dispositivos-medicos">https://www.infarmed.pt/web/infarmed/entidades/dispositivos-medicos/vigilancia-de-dispositivos-medicos</a>
Romania	<b>Oana Arsenescu Georgeta Herta</b> tel: +40 21.222.86.52, +40 21 260.01.58, +40 21 260.01.59 fax: +40 21.222.86.83 <b>National Agency for Medicines and Medical Devices of Romania</b> 58, Sos. Nicolae Titulescu, sector 1, Bucharest E-mail: mdevice@anm.ro, georgeta.herta@anm.ro	<a href="https://www.anm.ro/">https://www.anm.ro/</a>
Slovenija / Slovenia	<b>JAZMP - Agency for Medicinal Products and Medical Devices of the Republic of Slovenia</b> Slovenčeva ulica 22, SI - 1000 Ljubljana - - tel: +386 8 2000 500 E-mail: info@jazmp.si, meddev.vigilance@jazmp.si	<a href="https://www.jazmp.si/en/medical-devices/vigilance-of-medical-devices/">https://www.jazmp.si/en/medical-devices/vigilance-of-medical-devices/</a>
Slovenska Republika / Slovakia	<b>Darina Kaminská</b> - tel:+421 2 50701215, <b>State Institute for Drug Control, Medical Devices Section</b> Kvetna 11, SK - 825 08 Bratislava E-mail: darina.kaminska@sukl.sk	<a href="https://www.sukl.sk/hlavna-stranka/english-version?page_id=256">https://www.sukl.sk/hlavna-stranka/english-version?page_id=256</a>
Suomi / Finland	<b>Finnish Medicines Agency Fimea Medical Devices Unit</b> Mannerheimintie 166, P.O. box 55, FI-00034 FIMEA, FINLAND, tel:+358 29 522 3602 E-mail: meddev.vigilance@fimea.fi, registry@fimea.fi	<a href="https://www.valvira.fi/">https://www.valvira.fi/</a>
Sverige /	Swedish Medical Products Agency	<a href="https://www.lakemedelsverket.se/sv">https://www.lakemedelsverket.se/sv</a>

Country	Contact information	Website
Sweden	'Läkemedelsverket' Department of Medical Devices Box 26, SE-751 03 Uppsala, - tel:+46 18 174600, E-mail: meddev.central@lakemedelsverket.se -	
Norway / Norge	<b>AIMDD, MDD IVDMD - – Raymond Ludvigsen / Bjørn Kristian Berge Statens legemiddelverk/ Norwegian Medicines Agency</b> Postboks 6167 Etterstad, N-0602 Oslo, tel:+47 22 89 77 00 E-mail: meddev-no@legemiddelverket.no	<a href="https://legemiddelverket.no/English">https://legemiddelverket.no/English</a>
Iceland / Island	<b>Haukur Eggertsson</b> - tel:+354 520 2100, fax:+354 <b>561 2170 Icelandic Medicines Agency</b> Vínlandsleið 14, IS-113 Reykjavík E-mail: Haukur.Eggertsson@lyfjastofnun.is	<a href="https://www.lyfjastofnun.is/">https://www.lyfjastofnun.is/</a>
Liechtenstein	<b>Martin Stricker</b> , tel:+423 236 73 36 <b>Office of Public Health</b> Äulestrasse 51, Postfach 684, FL - 9490 Vaduz E-mail: medical.devices@llv.li -	<a href="https://www.llv.li/inhalt/1908/amtsstellen/amt-fur-gesundheit">https://www.llv.li/inhalt/1908/amtsstellen/amt-fur-gesundheit</a>
Turkey	<b>The Ministry Of Health, Turkish Medicine and Medical Device Agency</b> <b>Vice Presidency of Inspectorate, Medical Device Vigilance Unit</b> Söğütözü Mahallesi 2176. Sokak No:5 P.K. 06520 Çankaya Ankara/ TURKEY E-mail: meddev.vigilance@titck.gov.tr	/