Inogen One G5 / Inogen Rove 6

Column Change Guide - Multilingual Edition

Based on document 96-09240-00-01 D.

□□ English

INOGEN ONE G5 / Inogen Rove 6 - Column Change Guide

Fast and easy maintenance – the Inogen benefit.

1. Turn off the Inogen One G5/Rove 6 concentrator by pressing the power button.
2. Remove the concentrator from the carrying case, if attached.
3. Remove the battery from the concentrator.
4. Place the concentrator on its side so that the underside is visible. The metal column assembly can be seen on one side.
5. Unlock the column assembly by pushing the latch button away from the columns.
6. While holding the button open, slide the column assembly out by pulling on the column pull handle.
7. Remove the columns completely. Both columns are removed as one assembly.
8. If tape is applied to the column assembly, remove the tape.
9. Remove the dust caps from the new column assembly. Make sure there is no dust or debris where the dust caps were located.
10. Insert the column assembly into the concentrator. Do not leave the column ends exposed. Insert the assembly as soon as the dust caps have been removed.
11. Push the column assembly fully into the device. The spring-loaded latch button should fully return to the closed position.
12. Connect the AC power supply and plug it into an electrical outlet. Do not power on the concentrator.
13. Press and hold the + and – buttons for 5 seconds. The screen will display “sieve reset”. Release the buttons.
14. Press the alert/bell button once. The screen will display “sieve reset success.” Press the power button to turn on the concentrator and use normally.

NOTE / HINWEIS / REMARQUE / NOTA / HUOMAUTUS: Column change instructions are intended only when maintenance is required and are not intended for practice purposes.

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